







fluorochem.

New Product Launch & Offer

Thin Layer Chromatography (TLC) Silica Gel Plates



Aluminum Sheets

Why TLC?

- Easy to perform;
- Low cost;
- Small sample quantities needed;
- Possible recovery of the products;
- Short development time;
- Wide application.

ChemScene TLC Silica Gel Plates Advantages:

- Homogeneous, smooth and well adhering layers;
- Excellent separation efficiency;
- Consistent high quality from batch-to-batch and from plate-to-plate;
- Excellent customer service and product warranty.

Building Blocks / Pharmaceutical Intermediates / Chemical Reagents

CS-E0010 HPTLC Silica Gel Plate, Rock-Bottom Prices

Cat. No.	Length × Width	Layer Thickness	Qty/Pack
CS-E0009	2.5 × 7.5 cm	0.25 mm	320 Glass Plates
CS-E0010	2.5 × 5 cm	0.25 mm	640 Glass Plates
CS-E0011	5 × 20 cm	0.25 mm	80 Glass Plates
CS-E0013	20 × 20 cm	0.25 mm	20 Glass Plates
CS-E0007	20 × 20 cm	1 mm	20 Glass Plates
CS-E0008	20 × 20 cm	0.2 mm	25 Aluminum Sheets

30,000+

Building Blocks

5,000+

Chemical Reagents

5,000+

Pharmaceutical Intermediates